

## General Western Blot Protocol

1. Polyacrylamide gel electrophoresis and blotting
  - a. Add an appropriate amount of electrophoresis sample buffer (1X = 125mM Tris-HCl pH 6.8, 2% SDS, 5% glycerol, 0.003% bromophenol blue, and 1% beta-mercaptoethanol) to all samples.
  - b. Heat to 95° C for 3-5 minutes.
  - c. Load 5-100 ug total protein in a volume that is appropriate for the size of the wells.
  - d. Electrophorese proteins for the appropriate time according to the manufacturers' instructions.
  - e. Transfer proteins from the gel to a suitable membrane (e.g. nitrocellulose, PVDF, etc.) following the manufacturers protocol for transfer.
  - f. High molecular weight proteins (>120 kDa) can be wet transferred more efficiently if transfer time is increased and 0.05% SDS is included in the transfer buffer.
  
2. Blocking
  - a. Remove the filter from the transfer apparatus and rinse in PBST/TBST to remove loose acrylamide.
  - b. Transferred proteins can be visualized by staining the membrane for 15-30 minutes with Ponceau S.
  - c. Remove stain from filter by washing with PBST/TBST.
  - d. Place filter into blocking solution.
  - e. Block for 30 minutes at 37° C, 1 hour at room temperature, or overnight at 4° C.
  
3. Incubation with primary antibody
  - a. Decant the blocking buffer and add the antibody, diluted in blocking buffer as suggested in the product description sheet.
  - b. Incubate with agitation for 30 minutes at 37° C, 1 hour at room temperature, or overnight at 4° C.

4. Incubation with secondary antibody

- a. Wash for 30 minutes with agitation in wash buffer (TBS or PBS with 0.1% Tween 20), changing the wash buffer every 5 minutes.
- b. Decant the wash solution and add HRP-conjugated secondary antibody, diluted in 5% non-fat dry milk in wash buffer.
- c. Incubate for 30 minutes at 37° C, 1 hour at room temperature, or overnight at 4° C.
- d. Decant the antibody conjugate and wash for 30 minutes with agitation in wash buffer (TBS or PBS with 0.1% Tween 20), changing the wash buffer every 5 minutes.

5. Substrate incubation (ECL)

- a. Decant washing buffer and place the blot in a plastic bag or clean tray containing the development working solution (0.125 ml/cm<sup>2</sup>) for 1-5 minutes.
- b. Agitate bag or tray to cover the surface of the membrane.
- c. Remove the blot from the bag or tray and place it between two pieces of write-on transparency film.
- d. Smooth over the covered blot to remove air bubbles and excess substrate.
- e. Expose to X-ray film or any sensitive screen. (Check manufacturer's instructions for specific ECL reagents and procedures.)