

# NEUROMICS



## FluoGreen Tracer

## Protocol

**Catalog Number:** SF41000

**Applications:** This tracer can stain any animal cells including neurons, glia, lymphocytes and more.

*Store Frozen and Thaw Before Using.*

### Protocol

1. Calculate the volume of the working solution
2. Dilute FluoGreen Tracer stock solution with PBS (1:100)
3. Add diluted FluoGreen Tracer solution to cells or tissue sections and incubate from 1 to 5 min at room temperature
4. Rinse cells or tissue section with PBS
5. Mount using iBright Plus mounting media

| Name                         | Catalog #  | Type | Species        | Applications           | Size  | Price |
|------------------------------|------------|------|----------------|------------------------|-------|-------|
| <a href="#">i-BRITE Plus</a> | SF40000-1  | Kit  | Not Applicable | All Fluorescent Images | 1 ml  | \$29  |
| <a href="#">i-BRITE Plus</a> | SF40000-5  | Kit  | Not Applicable | All Fluorescent Images | 5 ml  | \$99  |
| <a href="#">i-BRITE Plus</a> | SF40000-10 | Kit  | Not Applicable | All Fluorescent Images | 10 ml | \$149 |

### Application Notes

Excitation and emission wavelengths of the fluorochrome:

Excitation is: 495 nm

Emission is: 510 nm

**FOR RESEARCH USE ONLY**

NEUROMICS' REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

**02/08v1**

[www.neuromics.com](http://www.neuromics.com)

Neuromics • 5325 West 74<sup>th</sup> Street, Suite 8 • Edina, MN 55439  
phone 866-350-1500 • fax 612-677-3976 • e-mail [pshuster@neuromics.com](mailto:pshuster@neuromics.com)