

NEUROMICS



E18 Primary Rat Neurospheres

Catalog Number: NS36100

Components : Approximately 1×10^6 cells

12 mls Culture Media: [Rat Neurosphere/Progenitor Proliferation Media \(Neuromics' Catalog#: M36111\).](#)

Cell stored in **2mls Hibernate EB:** Hibernate E®/B27®/Glutamax™
5 mls Hibernate E-Ca:

Description: Live, unseparated, fresh from E18 rat cortex/hippocampus including subventricular zone. One order is from one brain which will provide over 1 million neurospheres. NNSph can be used for expansion as neurospheres with repeated passage or differentiation into neurons or astrocytes.

Shipping/Storage: Neurospheres are shipped refrigerated. Cells are stable for up to 1 weeks when stored at 4-8°C. It is recommended to culture the cells as soon as possible after receiving cells according to Application Notes outlined below.

Application Notes

Materials Needed Not Provided

- Ultra Low Adhesion Wells (Corning)
- Papain (Worthington or Sigma P4762) for enzymatic dissociation
- Sterile pipette tips or sterile Pasteur pipette Sterile centrifuge tubes
- Centrifuge to operate at 200xg Water bath at 30°C Additional media:
 - [Rat Neurosphere/Progenitor Proliferation Media \(Neuromics' Catalog#: M36111\).](#)

Store tissue at 4°C until ready to use (note: use tissue right away for highest yield).

Preparations

1. Use Corning , 6-well ultra-low attachment plates (Fisher Scientific: 3471). No additional substrate is required.
2. Prepare cell dissociation solution by dissolving 6 mg sterile papain in 3 ml Hibernate E-Ca for a final working concentration of 2 mg/ml papain. Incubate for 10 min. at 30°C.
3. Fire polish the tip of a 9" salinized pasteur pipette to an opening of app. 0.5 mm.
4. Aliquot 20 µl Trypan Blue (Sigma T8154) into a 0.5 ml tube for step 9.

Dispersal (in a sterile hood)

1. With the salinized pasteur pipette, carefully transfer Hibernate EB (HEB) to a sterile tube (save for step 3) leaving the tissue with minimal HEB.
2. Add 2 ml of dissociation media to the tissue and incubate for 10 min. at 30°C. Gently swirl every 5 min.
3. Remove cell dissociation solution leaving the tissue at the bottom. Return the HEB medium from step 1.
4. With the salinized pasteur pipette, triturate tissue for about 1 min. (90% tissue dispersal) avoiding air bubbles.
5. Let the undispersed tissue settle for 1 min.
6. Transfer supernatant with cells to a new sterile 15 ml tube. Leave app. 5 µl of the HEB containing debris.
7. Spin 1100rpm (200 x G), 1 min. Discard supernatant leaving app. 5 µl of the HEB containing the pellet.
8. Disperse the pellet of cells (flick the bottom of the tube with a finger) and re-suspend the cells in 4 ml. of Neurosphere/Progenitor Proliferation Media
9. Aliquot 20 µl of cell solution into 0.5 ml. tube containing 20 µl Trypan Blue (1:2 dilution).
10. Count in hemacytometer: calculate cells/ml. Phase bright.

FOR RESEARCH USE ONLY V3-11/2012

Neuromics' reagents are for in vitro and certain non-human in vivo experimental use only and not intended for use in any human clinical investigation, diagnosis, prognosis, or treatment. We disclaim all liability in connection with the use of the information contained herein or otherwise, and all such risks are assumed by the user.

Cell Plating (at room temperature in a sterile hood)

1. Dilution cells in proliferation media and (2.0 ml/cm^2) and plate at $16,000 \text{ cells}/10 \text{ cm}^2$ or desired concentration.
2. Incubate 37°C , 5% CO_2 , 9% oxygen, 95 % humidity.
3. After 2-3 days neuroprogenitors are present.
4. Change $\frac{1}{2}$ the proliferation media every 3-4 days. Equilibrated at 37° 5% CO_2 .

Frequently Asked Questions**1. Why do I see clusters forming right away?**

Seeding cells at $30,000/\text{cm}^2$ will result in rapid cluster formation. Without attachment, they will stay as progenitors in Neuropro and grow into neurospheres that can be harvested in 4-7 days. Evaluation with immunostaining for nestin or other antibodies can begin after 4 days. To see clonal growth, you can plate the original sample or these neurospheres at limiting dilution of $50 \text{ cells}/\text{cm}^2$ of tissue culture plastic or ultra-low adhesion plastic in the same medium.

To see pluripotency, plate dissociated neurospheres at 5 to 15 thousand cells/ cm^2 on substrates coated with poly-d-lysine (50 $\mu\text{g}/\text{ml}$ water).

2. Should I seed single cells to do the neurosphere assay?

No, a neurosphere contains a cluster of cells. You can fix them and stain as a cluster or dissociate them with papain or trypsin, count to determine yield and/or fix to evaluate individual cells.

3. The cells were single after seeding but they formed cluster today. The sizes of some of the clusters are very big. Is this the right sign?

Yes.

4. Should I observe the neurosphere formation to evaluate the cells' progenitor stage?

No, wait at least 4 days.

5. How do I distinguish the neurosphere and cluster?

You can pass the neurospheres onto low adhesion plastic or onto adhesive surface and see differentiation. As long as cells are on low-adhesion plastic, they will stay as progenitors.

6. Do I need to do immunostaining at this stage?

Wait at least 4 days.

FOR RESEARCH USE ONLY V3-11/2012

Neuromics' reagents are for in vitro and certain non-human in vivo experimental use only and not intended for use in any human clinical investigation, diagnosis, prognosis, or treatment. We disclaim all liability in connection with the use of the information contained herein or otherwise, and all such risks are assumed by the user.