

Human NeuroProgenitors

Catalog #: NP4200, NP4200AA, NP4200OP
Storage: Liquid Nitrogen until ready to plate
Cell Count: 500,000 or 1,000,000

General Information

There are several challenges associated with In vitro neurological research. Specifically, there is a challenge in establishing batch-to-batch consistency, as well as high expansion capability. Induced pluripotent stem cell (iPSC)-derived neuronal progenitor cells have high yield capability and can differentiate into various neurological cells types, making them an ideal model for basic science, regenerative medicine and drug discovery. Our patient derived human neuronal progenitor cells were made from PBMCs, collected from substance abusing patients. The cells may be used for experimental applications including drug development, neurotoxicity, neurogenesis, and electrophysiology.




Image: Neuroprogenitors at 95% confluency

Characterization of the cells

Routine characterization of neuroprogenitors includes immunofluorescence staining and morphological observation of plated and suspension cells. Following adhesion to laminin and differentiation cells stain positive for Beta Tubulin III, specific for neuronal lineage, and Glial Fibrillary Acidic Protein (GFAP), denoting astrocyte lineage

All cells are performance assayed and test negative for HIV-1 and 2, mycoplasma, Hepatitis-A, B and C, bacteria, yeast and fungi. Cell viability, cell number and morphology are measured after recovery from cryopreservation.

Recommended Products

NeuroProgenitor Medium (cat# NM42400).

FOR RESEARCH USE ONLY

NEUROMICS REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. v1-09809

www.neuromics.com

Neuromics Antibodies • 5325 West 74th Street, Suite 8 • Edina, MN 55439
Phone 866-350-1500 • fax 612-677-3976 • email: pschuster@neuromics.com

Handling of Arriving Cells

Store in liquid nitrogen until ready to culture.

Note: Handling human derived products is potentially biohazardous. Always wear gloves and safety glasses when working these materials. Never mouth pipette. We recommend following the universal procedures for handling products of human origin as the minimum precaution against contamination.

Note: Always work under aseptic conditions.

CULTURE PROTOCOL

1. Prepare a Vitronectin (Gibco A14700) coated 35mm plate with a final concentration of 0.5ug/cm²-1.0ug/cm².
2. Rapidly thaw (<2 minute) frozen vial of NP-cells in a 37°C water bath. Remove vial from water bath just before the last trace of ice has melted. Spray vial with 70% ethanol to prevent contamination in tissue culture hood.
3. Thaw cells in 37°C water bath with agitation.
4. Wash once with 5 mLs NeuroProgenitor Medium
5. Pellet at 1000 rpm for 2 mins at room temp.
6. Resuspend cells in 2 mLs NeuroProgenitor and add to Vitronectin coated plate (created in step 1).
7. Place cells in 37°C incubator at 5% CO₂, replacing half of media every day until confluent (typically 2-3 days).
8. Cells should be passaged using Versene solution (Gibco). Wash cells with PBS (without Calcium and magnesium).
9. Add 1mL of Versene solution to cells (35mm plate). Place cells back into incubator for 2 mins.
10. Remove cells from incubator and place into tissue culture hood. Using a 1mL pipettor, pipette the 1mL solution (step 9) over the plate vigorously until cells are no longer adherent. Pipette solution, with cells, into a 15 mL centrifuge tube. It may be necessary to wash plate with additional PBS or media (NeuroProgenitor) in order to remove all cells. Adherent cell are visible to the eye, however it is recommended to examine the plate microscopically to ensure cells have been collected.
11. Pellet cells at 1000 rpm for 2 mins at room temperature.
12. Carefully remove supernatant without disturbing the cell pellet.
13. Resuspend cell pellet in 1 mL of NeuroProgenitor media and add half of the volume to two 35mm Vitronectine coated plates (previously coated).
14. Place cells in 37°C incubator at 5% CO₂, replacing half of media every day until confluent (typically 2-3 days).

Note: Passaging should be performed at 90% confluency or greater. If cells are allowed to remain at 100% confluency for multiple days may inhibit or greatly reduce cell growth after passaging. Cells should be passaged 1:2, and anything greater may be too sparse and cells may not recover, nor grow to confluence.

Note: Should any issues arise while using our cells, our team is here to help troubleshoot any issues. Our cells are backed by our one-time replacement or refund policy. Our recommended protocol including recommended products must be used to be eligible for replacement or refund. Cells that have been refrozen are no longer eligible for refund or replacement.

FOR RESEARCH USE ONLY

NEUROMICS REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. v1-09809

www.neuromics.com

Neuromics Antibodies • 5325 West 74th Street, Suite 8 • Edina, MN 55439
Phone 866-350-1500 • fax 612-677-3976 • email: pshuster@neuromics.com