


Human Cardiomyocytes– SC00A1-CM

Protocol for cell culture of human neural lineages derived from mesenchymal stem cells manufactured by Vitro Biopharma



Images: Left: Cardiomyocytes (SC00A1-CM- stained for cardiac tryponin T+ dapi). Right: Human Cardiomyocytes stained with Cardiac Tryponin T (green) and counter stained with DAPI (blue).

Establishing cultures from cryopreserved cells:

Use of these cell lines requires prior experience in standard methods of mammalian cell culture. In particular, sterile technique is required in a dedicated cell culture facility that is free from contamination.

Cryopreserved cells may be used to establish cultures immediately upon receipt or they may be stored for use at a later time. If stored, it is preferable to store in the vapor phase of liquid N₂. Storage in a –80°C freezer may be used but is likely to result in diminished cell viability proportional to storage time. NOTE: These cells can be cultured for 10 to 20 days and cryopreservation is not recommended. Do not thaw until ready for use.

To establish cultures from cryopreserved cells, first ensure that adequate equipment and reagents are available to perform the necessary procedures in a timely manner. These cells require culture in a 37°C, CO₂ cell culture incubator, calibrated to 5% CO₂ and 1% to 5% O₂. However, these cells can be cultured at ambient oxygen levels (~20% O₂) which results in reduced growth rates. Please contact technical services for information about various equipment and instrumentation options available to establish reduced oxygen cultures. Also, a water bath equilibrated to 37°C is needed. Required reagents include 1 x PBS, Fisher Catalog number BP665-1 or equivalent at room temperature or 37°C and Neural culture medium. We provide Cardiomyocyte maintaining media (PC00B6) for optimizing healthy cells and stability.

To establish cultures from frozen cell stock, it is first necessary to rapidly thaw cells at 37°C. Remove the desired number of vials containing cells from liquid nitrogen.

FOR RESEARCH USE ONLY

NEUROMICS' REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. V1-06/2012

www.neuromics.com

Neuromics Antibodies • 5325 West 74th Street, Suite 8 • Edina, MN 55439 phone 866-350-1500 • fax 612-677-3976 • e-mail pshuster@neuromics.com



Exposure of closed vials containing liquid nitrogen to a 37°C water bath is an explosion hazard! It is essential to ensure that no liquid nitrogen is present in the vial! If liquid nitrogen is present in the vial, allow this to evaporate before proceeding. Please use proper precautions including appropriate gloves to protect skin from exposure to liquid nitrogen, eye protection and other personal protective equipment when transferring vials containing cryopreserved cells from liquid nitrogen to a 37°C water bath.

Provide continuous agitation, e.g., swirling, to the vial while it is submerged in the 37°C water bath. Continue with agitation until the cells are completely thawed and no ice remains within the cell suspension, usually about 1 to 2 minutes. Maximum cell viability is dependent on rapid and complete thawing of frozen cells.

Count cells by a suitable method including a hemacytometer or automated cell counting device and determine the concentration of cells within the cell suspension. Our products are provided in 0.5 ml (500 l) of cryopreservation medium at a nominal concentration of 1×10^6 cells/ml. We recommend direct inoculation of cultures from the cell-cryopreservation media suspension. (Washout of the cryopreservative has been shown to decrease viability.)

We suggest establishing the culture at a plating density of about 50,000 to 75,000 cells/cm² in suitable tissue culture dishes or flasks. Please note that extracellular matrices are important in culture of all cells, especially attachment dependent cells. Differentiated neural lineages require laminin or fibronectin coated flasks. We recommend using Corning Biocoat Laminin Cellware (Catalog Number 354533). Add the appropriate volume of MSC culture medium to the plate or flask to be used for culture following the guidelines of the plate or flask manufacturer. We typically add 10 ml of medium per T-25. Use these guidelines to determine the appropriate volume of medium for your application. Following inoculation with the appropriate volume of cell suspension, gently agitate the flask or plate to ensure homogeneous distribution of the MSCs with the cell culture medium. Allow cultures to incubate in 5% CO₂ in ambient or reduced O₂ as noted above at 37°C in a humidified environment. Monitor cell growth by visual inspection. This should require about 10 to 20 days of continuous culture, but this time depends on several factors (e.g. incubation conditions, media used, Rock inhibitor use, etc). Thus the cultures should be monitored by inspection with an inverted microscope with appropriate magnification e.g., 100x. Once cells become attached, within 3-5 days, change the media every other day until desired confluency is achieved.

Note: Should any issues arise while using our cells, our team is here to help troubleshoot any issues. Our cells are backed by our one-time replacement or refund policy. Our recommended protocol including recommended products must be used to be eligible for replacement or refund. Cells that have been refrozen are no longer eligible for refund or replacement.

FOR RESEARCH USE ONLY

NEUROMICS' REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. V1-06/2012

www.neuromics.com

Neuromics Antibodies • 5325 West 74th Street, Suite 8 • Edina, MN 55439 phone 866-350-1500 • fax 612-677-3976 • e-mail pshuster@neuromics.com