

GFP-Immortalized Human Brain Microglia Cells

Catalog #: HBMCs001-GFP

Cell #: > 5X10⁵ cells

Storage: Liquid Nitrogen until ready for culture
While culturing keep in 37°C CO₂ Incubator

Product Format: Frozen vial

GENERAL INFORMATION

Human brain microglia cells are isolated from healthy human brain tissue. GFP-Immortalized Human Brain Microglia Cells are selected from HBMgs infected with lentiviruses expressing hTert and GFP with puromycin. HBMgs can be cultured in long-term (tested >20 passages). It is recommended to use Alpha-glia Expansion medium (AGEM-001) for the culturing of HBMgs.

Product is for Research use only.

Frozen vials are shipped in a dry ice package.

CHARACTERIZATION OF THE CELLS

- Cytoplasmic F4/80 Positive
- CD68 Positive
- HBMgs are tested negative for HIV-1, HBV, HCV, and mycoplasma

HANDLING OF ARRIVING CELLS

When you receive the cells in a frozen vial, you can transfer the vial of cells into a -80C freezer for short period storage or a liquid nitrogen tank for long-term storage. Thaw the cells in a 37C water bath, and then transfer the cells into a T25 flask pre-coated with Universal Coating Solution (AC002) as described in following details.

PROTOCOL FOR THAWING THE CELLS AND SUBCULTURE

Note: If you have any questions or need clarification regarding the protocol for culturing these cells, please reach out to Dr. Jensen Auguste at (978) 608-1766 with your questions before beginning.

1. Prepare a flask (T-25 flask is recommended) by adding 5 ml of Universal Coating Solution (AC002) to a T-25 flask, leave the flask in an incubator overnight (minimum one hour in 37C incubator).
2. Rinse the coated flask with sterile PBS twice and add 7 ml of Alpha-glia Expansion medium (AGEM-001) to the flask. Leave the flask in the hood and go to thaw the cells.
3. Warm Alpha-glia Expansion medium before thawing the cells.
4. Place the vial in a 37C water bath, hold and rotate the vial gently until the contents are completely thawed. Remove the vial from the water bath immediately, wipe it dry, rinse the vial with 70% ethanol and transfer it to a sterile field. Remove the cap, being careful not to touch the interior threads with fingers. Using 1 ml Eppendorf pipette gently resuspend the contents of the vial.
5. Dispense the contents of the vial into the neutralized coated culture vessels. A seeding density of more than 10,000 cells/cm² is recommended.
6. Replace the cap or cover, and gently rock the vessel to distribute the cells evenly. Loosen cap if necessary, to permit gas exchange.

FOR RESEARCH USE ONLY

NEUROMICS REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. v1-09809

www.neuromics.com

Neuromics Antibodies • 5325 West 74th Street, Suite 8 • Edina, MN 55439
Phone 866-350-1500 • fax 612-677-3976 • email: pshuster@neuromics.com

7. Return the culture vessels to the incubator (For best result, do not disturb the culture for at least 16 hours after the culture has been initiated. Change the growth medium the next day to remove the residual DMSO by centrifuging the unattached cells down, resuspending the cells with 5ml fresh medium, and add the cells back to the flask. Then change the medium every other day thereafter).
8. Change the medium every two to three days thereafter.

Note: Should any issues arise while using our cells, our team is here to help troubleshoot any issues. Our cells are backed by our one-time replacement or refund policy. Our recommended protocol including recommended products must be used to be eligible for replacement or refund. Cells that have been refrozen are no longer eligible for refund or replacement.

FOR RESEARCH USE ONLY

NEUROMICS REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. v1-09809

www.neuromics.com

Neuromics Antibodies • 5325 West 74th Street, Suite 8 • Edina, MN 55439
Phone 866-350-1500 • fax 612-677-3976 • email: pshuster@neuromics.com