



AlphaBioCoat Precoated 60 mm Dishes (Pack of 10)

Catalog #: AC001-60MM

Product Size: 60 mm Dish (Pack of 10)

Storage: 2-30°C

Shelf Life: 6 Months from Date of Receipt

GENERAL INFORMATION

AlphaBioCoat precoated 60 mm dishes have a uniform and consistent application of high-quality AlphaBioCoat solution (cat. AC001) on a clear polystyrene surface. Dishes have a raised ring around the circumference (SureGrip) of the dish that provides a secure grip of both the dish and lid. Embossed arrows on the dish and lid provide orientation for cell site observation. Dishes are aseptically processed and are non-pyrogenic. Dishes are available in quantities of 10 dishes per sleeve.

NOTE: Dishes need to be rinsed with 1X media or PBS prior to seeding cells.

Sterilization Method	Aseptically Processed
Stackable	Yes
Ventilation Cams	Yes
Flask Polymer	Polystyrene
Tissue Culture Treated Prior to Coating	Yes
Growth Area per Dish	22.9 cm ²
Typical Working Volume per Flask	5 ml

STORAGE AND USE

All procedures should be done under sterile conditions using aseptic techniques. Dishes must be rinsed with 1X media or PBS prior to seeding cells. Precoated flasks have a minimum shelf life of 6 months from the date of receipt and can be stored from 2-30°C. Flasks are shipped at room temperature.

FOR RESEARCH USE ONLY

NEUROMICS REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. v1-09809

www.neuromics.com

Neuromics Antibodies • 5325 West 74th Street, Suite 8 • Edina, MN 55439
phone 866-350-1500 • fax 612-677-3976 • e-mail pshuster@neuromics.com