

NEUROMICS

Caspase-12 Data Sheet

Catalog Number:	RA15048	Host:	Rabbit
Product Type:	Affinity purified	Species Reactivity:	Rat, Mouse
Immunogen Sequence:	Recombinant mouse Caspase-12, amino acids 1-419 with a carboxyl-terminal six histidine tag.	Format:	Liquid 1mg/ml 0.2 µm filtered solution in phosphate-buffered saline (PBS) with 5% trehalose.
Applications:	Western Blot: 0.25 µg/mL Tested in Western blot using mouse CTLL-2 and C2C12 cells and rat PC12 cells. Dilutions listed as a recommendation. Optimal dilution should be determined by investigator.		
Storage:	Antibody can also be aliquotted and stored frozen at -20° C to -70° C in a manual defrost freezer for six months without detectable loss of activity. The antibody can be stored at 2° - 8° C for 1 month without detectable loss of activity. <i>Avoid repeated freeze-thaw cycles.</i>		

Application Notes

Specificity

The antibody detects endogenous rat and mouse Caspase-12 proenzyme.

Western blotting

Blotting Buffer	Blocking Solution	Antibody Solution
25 mM Tris, pH 7.4 0.15 M NaCl 0.1% Tween 20	5% nonfat dry milk in Blotting Buffer, Adjust pH to 7.4	5% nonfat dry milk in Blotting Buffer Adjust pH to 7.4

1. Transfer the electrophoresed proteins to Immobilon membrane (Millipore) and incubate the membrane for 1 hour at room temperature in Blocking Solution.
2. Incubate the membrane overnight at 4° C in Antibody Solution containing 0.25 µg/mL anti-Caspase-12.
3. Wash the membrane at room temperature for 1 hour with 5 or more changes of Blotting Buffer. Changing the membrane containers often reduces background.
4. Incubate the membrane at room temperature for 1 hour in Antibody Solution containing a HRP-conjugated anti-rabbit IgG secondary antibody.
5. Wash the membrane for 1 hour with 5 or more changes of Blotting Buffer.
6. Detect with ECL Reagent (Amersham).

Cell lysates for Western blots: To prepare total cell lysates, cells are solubilized in hot 2X SDS gel sample buffer (20 mM dithiothreitol, 6% SDS, 0.25 M Tris, pH 6.8, 10% glycerol, 10 mM NaF and bromophenyl blue) at 2×10^6 - 1×10^7 cells per mL. The extracts are heated in a boiling water bath for 5 minutes and then sonicated with a probe sonicator with 3 - 4 bursts of 5 - 10 seconds each. Samples are diluted with 1X SDS sample buffer to the desired concentration.

FOR RESEARCH USE ONLY

NEUROMICS' REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

www.neuromics.com

Neuromics Antibodies • 11200 Hampshire Avenue South • Bloomington, MN 55438
phone 507-645-8020 • fax 612-677-3976 • e-mail pshuster@neuromics1.com

10/04v1