



Catalog Number:	RA15043	Host:	Rabbit
Product Type:	KLH coupled synthetic peptide CADLREDPDRQDHHPGSGAQ. The sequence corresponds to the carboxyl terminus of Ubiquitin+1 (van Leeuwen, F.W. et al., 1998, Science 279:242). Cysteine was added to the amino-terminus for conjugation to KLH and for coupling to an affinity matrix.	Species Reactivity:	Human
Immunogen Sequence:	Hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with rat brain synaptic plasma membranes.	Format:	Liquid 1mg/ml Solution in phosphate-buffered saline (PBS) with 5% Trehlose
Applications:	Immunohistochemistry -3-25 µg/mL Immunoprecipitation-See reference Western Blot-See reference		
	Dilutions listed as a recommendation. Optimal dilution should be determined by investigator.		
Storage:	Antibody can be aliquotted and stored frozen at -20° C to -70° C in a manual defrost freezer for six months without detectable loss of activity. The antibody can be stored at 2° - 8° C for 1 month without detectable loss of activity. <i>Avoid repeated freeze-thaw cycles.</i>		

Application Notes

Specificity

The antibody detects human Ubiquitin+1. It does not cross-react with Ubiquitin.

Western blotting

Blotting buffer
25 mM Tris, pH 7.5
0.15 M NaCl
0.05% Tween 20

Blocking solution
2% nonfat dry milk in
blotting buffer
pH to 7.5

Antibody solution
1% nonfat dry milk in
blotting buffer
pH to 7.5

1. Transfer the electrophoresed proteins to Immobilon filters (Millipore) and incubate the membrane for 1 hour at room temperature in blocking solution.
2. Incubate the membrane overnight in antibody solution containing 1.0 µg/mL rabbit anti-Ubiquitin+1.
3. Wash the membrane at room temperature for 1 hour with 5 or more changes of blotting buffer. Changing the membrane containers often reduces background.
4. Incubate the membrane for 1 hour at room temperature in antibody solution containing a 1:2,000 dilution of HRP-conjugated Protein A (Amersham).
5. Wash the membrane for 1 hour with 5 or more changes of blotting buffer.
6. Detection was with ECL Reagent (Amersham).

FOR RESEARCH USE ONLY

NEUROMICS' REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

11/10v2

www.neuromics.com

Neuromics • 5325 West 74th Street, Suite 8 • Edina, MN 55439
phone 866-350-1500 • fax 612-677-3976 • e-mail pshuster@neuromics1.com

Cell lysates for western blottings: To prepare total cell lysates, cells are solubilized in hot 2x SDS gel sample buffer (20 mM dithiothreitol, 6% SDS, 0.25 M Tris, pH 6.8, 10% glycerol, and bromophenyl blue) at 2×10^6 - 1×10^7 cells per mL. The extracts are heated in a boiling water bath for 5 minutes and then sonicated with a probe sonicator with 3 - 4 bursts of 5 - 10 seconds




Image: Immunoblots of 2 ng recombinant human Ubiquitin and 2, 1, and 0.5 ng recombinant human Ubiquitin+1. Samples were electrophoresed on 15% gels and immunoblotting was with 1.0 $\mu\text{g}/\text{mL}$ anti-Ubiquitin+1. A one minute exposure is shown.

FOR RESEARCH USE ONLY

NEUROMICS' REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RSKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

11/10v2

www.neuromics.com

Neuromics • 5325 West 74th Street, Suite 8 • Edina, MN 55439
phone 866-350-1500 • fax 612-677-3976 • e-mail pshuster@neuromics1.com