



Catalog Number:	RA15050	Host:	Goat
Product Type:	Affinity Purified Antibody	Species Reactivity:	Human
Immunogen Sequence:	KLH coupled synthetic peptide, MFQIPEFEPSEQEDSSSAERGC (corresponding to amino acids 1 - 21 of human Bad). Cysteine was added to the carboxyl-terminal for coupling to an affinity matrix.	Format:	Liquid 1mg/ml Solution in phosphate-buffered saline (PBS) with 5% Trehlose
Applications:	Western Blot-1 µg/mL		
	Dilutions listed as a recommendation. Optimal dilution should be determined by investigator.		
Storage:	Antibody can be aliquotted and stored frozen at -20° C to -70° C in a manual defrost freezer for six months without detectable loss of activity. The antibody can be stored at 2° - 8° C for 1 month without detectable loss of activity. <i>Avoid repeated freeze-thaw cycles.</i>		

Application Notes

Specificity

The antibody is known to react with human and mouse Bad.

Western Blot

An antibody concentration of 1.0 µg/mL is recommended.

Western Blot Protocol

Blotting Buffer

25 mM Tris, pH 7.5
0.15 M NaCl
0.05% Tween 20

Blocking Solution

5% nonfat dry milk in blotting buffer
pH 7.5

Antibody Solution


2% nonfat dry milk in blotting buffer
pH 7.5

1. Transfer the electrophoresed proteins to Immobilon filters (Millipore) and incubate the membrane for 1 hour at room temperature in blocking solution.
2. Incubate the membrane overnight at 4° C in antibody buffer containing 1.0 µg/mL rabbit anti-Bad.
3. Wash the membrane at room temperature for 1 hour with 5 or more changes of blotting buffer.
4. Incubate the membrane for 1 hour at room temperature in antibody buffer containing a 1:2,000 dilution of HRP-conjugated Protein A (Amersham).
5. Wash the membrane for 1 hour with 5 or more changes of blotting buffer.
6. Detection was with ECL Reagent (Amersham).

FOR RESEARCH USE ONLY

NEUROMICS' REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RSKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. V3-01/13

Cell lysates for Western blottings: To prepare total cell lysates, cells are solubilized in hot 2x SDS gel sample buffer (20 mM dithiothreitol, 6% SDS, 0.25 M Tris, pH 6.8, 10% glycerol, and bromophenyl blue) at 2×10^6 - 1×10^7 cells per mL. The extracts are heated in a boiling water bath for 5 minutes and then sonicated with a probe sonicator with 3 - 4 bursts of 5 - 10 seconds each. Samples are diluted with 1x SDS sample buffer to the desired concentration.



Immunoblots of SDS-extracts from human Jurkat (3×10^5), A431 (2×10^5), mouse CTLL2 (2.5×10^5), and L929 (2×10^5) cells. Extracts were blotted with 1.0 and 2.0 $\mu\text{g}/\text{mL}$ anti-Bad. Incubation with anti-Bad was overnight at 4°C and detection was by the ECL procedure (Amersham). One and five minute exposures are shown. The band that co-migrates with recombinant Bad is indicated by the arrow. The identity of the other bands is unknown.

FOR RESEARCH USE ONLY

NEUROMICS' REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RSKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. V3-01/13

www.neuromics.com

Neuromics • 5325 West 74th Street, Suite 8 • Edina, MN 55439
 phone 866-350-1500 • fax 612-677-3976 • e-mail pshuster@neuromics.com