



AlphaBioCoat Precoated 48-Well Plates (Pack of 10)

Catalog #: AC001-48

Product Size: 48-Multiwell Plate (Pack of 10)

Storage: 2-30°C

Shelf Life: 6 Months from Date of Receipt

GENERAL INFORMATION

AlphaBioCoat precoated 48-Well Plates with a flat bottom have a uniform and consistent application of high-quality AlphaBioCoat solution (cat. AC001) on a clear polystyrene surface. Plates are compatible for use with all common dispensers, readers, and washers. Plates are aseptically processed and are non-pyrogenic. Plates are available in quantities of 10 plates per sleeve.

NOTE: Plates need to be rinsed with 1X media or PBS prior to seeding cells.

| | |
|---|-----------------------|
| Sterilization Method | Aseptically Processed |
| Flask Polymer | Polystyrene |
| Tissue Culture Treated Prior to Coating | Yes |
| Growth Area per Well | 0.76 cm ² |
| Nominal Volume per Well | 1.6 ml |
| Typical Working Volume per Well | 0.3 ml |

STORAGE AND USE

All procedures should be done under sterile conditions using aseptic techniques. Plates must be rinsed with 1X media or PBS prior to seeding cells. Precoated plates have a minimum shelf life of 6 months from the date of receipt and can be stored from 2-30°C. Plates are shipped at room temperature.

FOR RESEARCH USE ONLY

NEUROMICS REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. v1-09809

www.neuromics.com

Neuromics Antibodies • 5325 West 74th Street, Suite 8 • Edina, MN 55439
phone 866-350-1500 • fax 612-677-3976 • e-mail pshuster@neuromics.com