



Epithelial Cryopreserve Media

Catalog #: FM301

Size: 20 ml

Storage: Store the medium at 2-8°C short term and -20°C long term

Format: Liquid

GENERAL INFORMATION

Epithelial Cryopreserve Media is qualified for use with human epithelial cells, it is optimized to allow for consistent cryopreservation and high viability upon thawing and plating.

Product is for Research use only. Our products are not authorized for human use, in vitro diagnostic procedures, or for therapeutic procedures.

Storage Condition: Store the medium at 2-8°C short terms and -20°C long term. Item is shipped with gel packs.

Note: To ensure sterility after 2 weeks or if there is concern that sterility was compromised during the supplementation process, the prepared medium may be re filtered with a 0.2 um filter.

CYROPRESERVATION PROTOCOL

1. Thaw cell culture freezing medium completely and mix well by gently swirling bottle. Keep freezing medium on ice during use.
2. Cells to be frozen should be in late log phase growth.
3. Monolayers will need to be dissociated. After dissociation, cells are resuspended in Human Epithelial Cell Basal Medium and counted to determine viability and number.
4. Centrifuge cells at 300 xg for 3 min. Remove the medium above the pellet.
5. Resuspend the cells in Epithelial Cryopreserve Media at a concentration of $\sim 4 \times 10^6$ cells/ml. Freeze 1 ml of cells/vial. After the cells have been resuspended and aliquoted into appropriate cryogenic storage vials, they can be placed in a freezing container and the normal freeze down procedure should be started within five minutes.
6. Cells must be stored at or below -80°C. For long term storage the cells should be stored in ultra-low temperature freezer (-150°C) or in liquid nitrogen (-196°C).

FOR RESEARCH USE ONLY

NEUROMICS REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. v1-09809

www.neuromics.com

Neuromics Antibodies • 5325 West 74th Street, Suite 8 • Edina, MN 55439
phone 866-350-1500 • fax 612-677-3976 • e-mail pshuster@neuromics.com