

# Survivin

# Data Sheet

**Catalog Number:** RA19022

**Product Type:** Rabbit IgG

**Immunogen Sequence:** EETAKKVRRAIEQLAAMD  
(a.a. 125-142) of human  
Survivin

**Format:** Affinity Purified  
liquid

**Species Reactivity:** Human, Mouse, Rat

**Applications:** Western blotting

**Dilutions:** WB 1:500 - 1:1000

**Storage:** Maintain at +2-8 degrees C for 3 months or at -20 degrees C for longer periods. Stable for 1 year. Avoid repeated freeze-thaw cycles.

## Application Note

Anti-Survivin antibody is at a concentration of 1 mg/ml in PBS containing 0.02% sodium azide.