

NEUROMICS



Goat anti-Chicken IgY–Unlabeled Data Sheet

Catalog Number:	CH23100	Host:	Goat
Product Type:	Affinity Purified	Species	Chicken
Description:	Affinity purified goat anti-chicken IgY, 1000 ul, 1.0 mg/ml-unlabeled	Reactivity:	
Applications:	Immunocytochemistry: 1:250 Immunohistochemistry: 1:250 Western blot: 1:5000 Dilutions listed as a recommendation. Optimal dilution should be determined by investigator.	Format:	Liquid PBS (pH 7.2; 10 mM; isotonic 0.9%, w/v) with sodium azide (0.02%, w/v).
Storage:	Store at 4°C in the dark. Under these conditions, the antibody should have a shelf life of at least 12 months (provided they remain sterile). Do not freeze the antibody unless you want to store it for longer periods of time. Note, however, that each time an antibody preparation is frozen, about half its binding activity is lost.		

Application Notes

FOR RESEARCH USE ONLY

Neuromics' reagents are for in vitro and certain non-human in vivo experimental use only and not intended for use in any human clinical investigation, diagnosis, prognosis, or treatment. We disclaim all liability in connection with the use of the information contained herein or otherwise, and all such risks are assumed by the user.

v1