



Catalog Number: CH23017

Product Type: Affinity Purified

Immunogen Sequence: Peptide corresponded to a region of the Coronin 1A gene sequence shared between the mouse(NM_009898) and human (NM_007074) gene products.

Host: Chicken

Species Reactivity: Human, Mouse

Format: 100 ul at a concentration Of 100 ug/ml in phosphate-buffered isotonic saline (PBS) with 0.02% sodium azide as an anti-microbial.

Applications: **Immunocytochemistry:** 1:1000
Immunohistochemistry: 1:1000

Dilutions listed as a recommendation. Optimal dilution should be determined by investigator.

Storage: **Store at 4°C in the dark.** Under these conditions, the antibody should have a shelf life of at least 12 months (provided they remain sterile). Do not freeze the antibody unless you want to store it for longer periods of time. Note, however, that each time an antibody preparation is frozen, about half its binding activity is lost.

Application Notes