



Recombinant Human IL-17

Datasheet

Catalog Number: PR27113

Product Type: Recombinant Protein

Source: *E.Coli*

Amino Acid Sequence: The sequence of the first five N-terminal amino acids was determined and was found to be, Gly-Ile-Thr-Ile-Pro.

Description/Molecular Mass: Interleukin-17A Human Recombinant produced in E.Coli is a homodimeric, non-glycosylated polypeptide chain containing a total of 264 amino acids and having a molecular mass of 30,250 Dalton. The IL-17 is purified by proprietary chromatographic techniques.

Protein Content: Protein quantitation was carried out by two independent methods:

1. UV spectroscopy at 280 nm using the absorbency value of 1.13 as the extinction coefficient for a 0.1% (1mg/ml) solution. This value is calculated by the PC GENE computer analysis program of protein sequences (IntelliGenetics).

2. Analysis by RP-HPLC, using a standard solution of IL-17 as a Reference Standard.

Purity: Greater than 98.0% as determined by:

- (a) Analysis by RP-HPLC.
- (b) Analysis by SDS-PAGE.

Biological Activity: The ED50 as determined by the dose-dependant induction of IL-6 in primary human foreskin fibroblasts was found to be 2 ng/mL, corresponding to a specific activity of 5×10^5 Units/mg.

Format: Lyophilized from a concentrated (1mg/ml) solution containing no additives.

Reconstitution: It is recommended to reconstitute the lyophilized Interleukin-17 in sterile 18MΩ-cm H₂O not less than 100µg/ml, which can then be further diluted to other aqueous solutions.

Storage: Lyophilized Interleukin-17 although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution IL-17 should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA).

FOR RESEARCH USE ONLY

NEUROMICS' REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

03/08v1

www.neuromics.com

Neuromics • 5325 West 74th Street, Suite 8 • Edina, MN 55439
phone 866-350-1500 • fax 612-677-3976 • e-mail pshuster@neuromics1.com