

NEUROMICS

MMP-24 (MT5-MMP) Data Sheet

Catalog Number:	RA15045	Host:	Rabbit
Product Type:	Protein A purified	Species Reactivity:	Human
Immunogen Sequence:	<i>E. coli</i> -derived GST-human MMP-24 (Tyr 125 to Cys 538) fusion protein	Format:	Liquid 1mg/ml 0.2 µm filtered solution in phosphate-buffered saline (PBS).
Applications:	Immunohistochemistry: 15 µg/mL Western Blot: 1 µg/mL Dilutions listed as a recommendation. Optimal dilution should be determined by investigator.		
Storage:	Antibody can also be aliquotted and stored frozen at -20° C to -70° C in a manual defrost freezer for six months without detectable loss of activity. The antibody can be stored at 2° - 8° C for 1 month without detectable loss of activity. <i>Avoid repeated freeze-thaw cycles.</i>		

Application Notes

Immunohistochemistry

For immunohistochemistry, the antibody will detect human MMP-24 in paraffin-embedded tissue sections. The working dilution is 15 µg/mL

FOR RESEARCH USE ONLY

NEUROMICS' REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

www.neuromics.com

Neuromics Antibodies • 210 Orchard St N • Northfield, MN 55057
phone 507-645-8020 • fax 612-677-3976 • e-mail psmuster@neuromics1.com