

Neuronal Nitric Oxide

Data Sheet

Catalog Number: RA19011

Product Type: Rabbit IgG

Immunogen Sequence: PRKGPPLPNGDTEC (a.a. 19-31) of neuronal NOS

Format: Affinity Purified liquid

Species Reactivity: Human

Applications: Western blotting

Dilutions: WB 1:500 - 1:1000

Storage: Maintain at +2-8 degrees C for 3 months or at -20 degrees C for longer periods. Stable for 1 year. Avoid repeated freeze-thaw cycles.
