

NEUROMICS



Recombinant Cytomegalo Virus Mosaic

Datasheet

Catalog Number: PR27123

Product Type: Recombinant Protein

Source: *E.Coli*

Description: The *E.Coli* derived recombinant artificial mosaic protein contains the CMV gB immunodominant regions, 60-81aa of gB protein strain AD169 repeated three times.

CMV Mosaic was purified by proprietary chromatographic technique.

Purity: CMV Mosaic protein is >95% pure as determined by 10% PAGE (coomassie staining).

Specificity: Immunoreactive with sera of CMV-infected individuals.

Applications: CMV Mosaic antigen is suitable for ELISA and Western blots, excellent antigen for detection of CMV with minimal specificity problems.

Format: 25mM Tris-Hcl pH 8.0, 1mM EDTA and 50% glycerol.

Storage: CMV Mosaic Protein is shipped at ambient temperature. Upon arrival, Store at -20°C

Stability is guaranteed for five years if frozen and one month in solution at room temperature.

FOR RESEARCH USE ONLY

NEUROMICS' REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RSKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

05/08v1

www.neuromics.com

Neuromics • 5325 West 74th Street, Suite 8 • Edina, MN 55439
phone 866-350-1500 • fax 612-677-3976 • e-mail pshuster@neuromics1.com