



Kidins220 (ARMS)

Data Sheet

Catalog Number:	RA19019	Host:	Rabbit
Product Type:	Affinity Purified Antibody	Species Reactivity:	Human, Mouse
Immunogen Sequence:	Reacts with residues 1747-1762[ASSESTGFGEERESIL] of the human 220kDa Kidins220 protein.	Format:	Liquid 1mg/ml in PBS containing 0.02% sodium azide.
Applications:	Western blot: 1:500-1:1000 Dilutions listed as a recommendation. Optimal dilution should be determined by investigator.		
Storage:	Maintain at +2-8 °C for 3 months or at -20 °C for longer periods. Stable for 1 year. Serum may have become trapped in top of vial during shipping. Centrifugation of vial is recommended before opening. Stable for at least 6 months at -20°C. Repeated freeze/thaw cycles compromise the integrity of the antiserum.		
References:	Iglesias T et al. J Biol Chem 2000 Dec 22;275(51):40048-56.		

Application Notes

FOR RESEARCH USE ONLY

Neuromics' reagents are for in vitro and certain non-human in vivo experimental use only and not intended for use in any human clinical investigation, diagnosis, prognosis, or treatment. We disclaim all liability in connection with the use of the information contained herein or otherwise, and all such risks are assumed by the user.

v1