

NEUROMICS



Neuropeptide Q2/NPQ2

Data Sheet

Catalog Number:	GP14109	Host:	Guinea Pig
Product Type:	Whole Serum	Species	Rat
Immunogen Sequence:	[NH ₂ -APQRLSE]-[K]7-bAla-OH	Reactivity:	
		Format:	100 ul. in PBS containing 0.05% sodium azide.
Applications:	Immunohistochemistry: 1:100-600 Dilutions listed as a recommendation. Optimal dilution should be determined by investigator.		
Storage:	Maintain at +2-8°C for 3 months or at -20°C for longer periods. Stable for 1 year. <i>Avoid repeated freeze-thaw cycles.</i>		

Application Notes

To improve staining intensity, it may be necessary to employ HRP-DAB (Vector Elite ABC kit) detection or enhance fluorescence signal using Tyramide Signal Amplification procedure.

FOR RESEARCH USE ONLY

v1

Neuromics' reagents are for in vitro and certain non-human in vivo experimental use only and not intended for use in any human clinical investigation, diagnosis, prognosis, or treatment. We disclaim all liability in connection with the use of the information contained herein or otherwise, and all such risks are assumed by the user.

phone 866-350-1500 • fax 612-677-3976 • www.neuromics.com