



Tyrosine Hydroxylase

Data Sheet

Catalog Number:	CH23006	Host:	Chicken
Product Type:	Affinity Purified	Species Reactivity:	Human, Mouse
Immunogen Sequence:	Peptide corresponded to a region of the Tyrosine Hydroxylase gene product shared between mouse (P24529, NCBI) and human (P07101, NCBI) gene products.	Format:	Liquid PBS (pH 7.2; 10 mM; isotonic 0.9%, w/v) with sodium azide (0.02%, w/v).
Applications:	Immunocytochemistry: 1:100 -1:200 Immunohistochemistry: 1:100 -1:200 Western blot: 1:250 -1:500 Dilutions listed as a recommendation. Optimal dilution should be determined by investigator.		
Storage:	Store at 4°C in the dark. Under these conditions, the antibody should have a shelf life of at least 12 months (provided they remain sterile). Do not freeze the antibody unless you want to store it for longer periods of time. Note, however, that each time an antibody preparation is frozen, about half its binding activity is lost.		

Application Notes

FOR RESEARCH USE ONLY

Neuromics' reagents are for in vitro and certain non-human in vivo experimental use only and not intended for use in any human clinical investigation, diagnosis, prognosis, or treatment. We disclaim all liability in connection with the use of the information contained herein or otherwise, and all such risks are assumed by the user.

v1