

# NEUROMICS



## HAP40 (Huntingtin-associated protein 40) Data Sheet

---

**Catalog Number:** RA19010

**Host:** Rabbit

**Product Type:** Affinity Purified Antibody

**Species Reactivity:** Mouse

**Immunogen Sequence:** DGHGQDTSGQLP (a.a. 314-325 ) of the 40 kDa mouse HAP40 protein

**Format:** Liquid. 100 ug in 100 ul (1 mg/ml) in PBS containing 0.02% sodium azide

**Applications:** Western blot: 1:500

Dilutions listed as a recommendation. Optimal dilution should be determined by investigator.

**Storage:** Maintain at +2-8°C for 3 months or at -20°C for longer periods. Stable for 1 year. *Avoid repeated freeze-thaw cycles.*

---

### Application Notes

#### *FOR RESEARCH USE ONLY*

NEUROMICS REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. v3-01/13

---

**www.neuromics.com**

Neuromics Antibodies • 5325 West 74<sup>th</sup> Street, Suite 8 • Edina, MN 55439  
phone 866-350-1500 • fax 612-677-3976 • e-mail [pshuster@neuromics.com](mailto:pshuster@neuromics.com)