

# CLAC-P (NC2-1 region)

# Data Sheet

---

**Catalog Number:** RA19007

**Product Type:** Rabbit IgG

**Immunogen Sequence:** NHGFLSADQQLIK (a.a. 171-183); NC2-1 region of human and mouse CLAC-P

**Format:** Affinity Purified liquid

**Species Reactivity:** Human

**Applications:** Western blotting

**Dilutions:** WB 1:500

**Storage:** Maintain at +2-8 degrees C for 3 months or at -20 degrees C for longer periods. Stable for 1 year. Avoid repeated freeze-thaw cycles.

---