



Catalog Number:	CH23101	Host:	Goat
Product Type:	Affinity Purified	Species:	Chicken
Description:	Affinity purified fluorescein-labeled goat anti-chicken IgY, 500 ul, 1.0 mg/ml. Antibody type: Polyclonal	Reactivity:	
Applications:	Immunocytochemistry: 1:250 Dilutions listed as a recommendation. Optimal dilution should be determined by investigator.	Format:	Liquid PBS (pH 7.2; 10 mM; isotonic 0.9%, w/v) with sodium azide (0.02%, w/v).
Storage:	Store at 4°C in the dark. Under these conditions, the antibody should have a shelf life of at least 12 months (provided they remain sterile). Do not freeze the antibody unless you want to store it for longer periods of time. Note, however, that each time an antibody preparation is frozen, about half its binding activity is lost.		

Application Notes