



Catalog Number:	MO25017	Host:	Mouse
Product Type:	Monoclonal Antibody	Species Reactivity:	Human and all mammalian
Immunogen Sequence:	IgG1 Native Uchl1 (PGP9.5) protein from brain was used as immunogen to generate the antibody (Day & Thompson, 1986).	Format:	500 ul. Tissue culture supernatant
Applications:	Western Blot 1:10000 Immunohistochemistry 1:5000 Immunohistochemistry-Paraffin 1:5000		


Dilutions listed only as a recommendation. Optimal dilution should be determined by investigator.

Storage: Store frozen. Aliquot as undiluted antisera and immediately place at -20°C. Antisera may have become trapped in top of vial during shipping. Centrifugation of vial is recommended before opening. Stable for at least 6 months at -20°C. Repeated freeze/thaw cycles compromise the integrity of the antiserum.

Application Notes

Description/Data:

Image: Mixed cultures of rat CNS cells stained with alpha-Internexin (red) and [chicken antibody to Microtubule associated protein 2](#) (MAP2- green). The a-interneixin is localized primarily in neuronal axons in these cultures, while the perikarya and dendrites of neurons stain strongly for MAP2.



FOR RESEARCH USE ONLY

NEUROMICS REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

www.neuromics.com

Neuromics Antibodies • 5325 West 74th Street, Suite 8 • Edina, MN 55439
phone 866-350-1500 • fax 612-677-3976 • pshuster@neuromics.com