



---

<b>Catalog Number:</b>	RA34005	<b>Host:</b>	Rabbit
<b>Product Type:</b>	Affinity Purified Polyclonal Antibody	<b>Species Reactivity:</b>	Human, Mouse, Rat
<b>Immunogen Sequence:</b>	Synthesized peptide derived from human Tau around the phosphorylation site of T534. Immunogen Region:150-230aa Gene ID: 4137; Swiss-Prot Key: P10636	<b>Format:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 100 ul aliquot Concentration: 1 mg/ml
<b>Applications:</b>	Phospho-Tau (T534) Polyclonal Antibody detects endogenous levels of Tau protein only when phosphorylated at T534. Western Blot: 1:500-2,000 Immunohistochemistry: 1:100-300 ELISA: 1:5,000		
<b>Storage:</b>	Dilutions listed as a recommendation. Optimal dilution should be determined by investigator. Antibody can also be aliquoted and stored frozen at -20° C in a manual defrost freezer for six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.		

---

### Application Notes

#### Description/Data:

Specificity: Phospho-Tau (T534) Polyclonal Antibody detects endogenous levels of Tau protein only when phosphorylated at T534.

*FOR RESEARCH USE ONLY*

NEUROMICS' REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RSKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.-V2/08/2012