

# Neuropeptide FF1

# Data Sheet

---

**Catalog Number:** RA19017

**Product Type:** Rabbit IgG

**Immunogen Sequence:** N terminal sequence  
MEGEP SQPPNSSWPLS  
and the C terminal sequence  
CSHLPLTIPAWDI of the  
human NPFF1 protein

**Format:** Affinity Purified  
liquid

**Species Reactivity:** Human

**Applications:** Western blotting

**Dilutions:** WB 1:500 - 1:1000

**Storage:** Maintain at +2-8 C for 3 months or at -20 C for longer periods. Stable for 1 year. Avoid repeated freeze-thaw cycles

---