

NEUROMICS



ENDO-Growth Media

Catalog #: MED001

Size: 500 ml

Storage: 2°C-8°C for up to 2 months
-20°C long term storage

Product Format: Liquid

General Information

ENDO-Growth Media contains two components freshly mixed for your convenience: 475 ml of ENDO-Basal media and 25 ml of ENDO-Growth supplement.

ENDO-Basal Media

Media containing central vitamins and amino acids, and 0.5% FBS.

ENDO-Growth Supplement

Contains the essential nutrients for endothelial cell proliferation; including 5% fetal bovine serum, recombinant growth factors and 1x penicillin and streptomycin

Shipping and Storage

Product is shipped with ice packs. It can be stored at 2°C-8°C for up two months and for long term storage keep at -20°C.

Protocol

1. Take the desired amount of ENDO-Growth Media need and warm to 37°C in a water bath or incubator
2. Avoid frequent temperature changes to the entire bottle of medium
3. The ENDO-Growth Media is stable for at least one month stored at 4°C
4. It is recommended that media is changed every day for regular cell culture

FOR RESEARCH USE ONLY

NEUROMICS REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. v1-09809

www.neuromics.com

Neuromics Antibodies • 5325 West 74th Street, Suite 8 • Edina, MN 55439
Phone 866-350-1500 • fax 612-677-3976 • email: pshuster@neuromics.com