



Peripherin

Data Sheet

Catalog Number:	CH23016	Host:	Chicken
Product Type:	Affinity Purified	Species:	Human, Mouse
Immunogen Sequence:	Recombinant protein expressed in bacteria	Reactivity:	
		Format:	Liquid PBS (pH 7.2; 10 mM; isotonic 0.9%, w/v) with sodium azide (0.02%, w/v).
Applications:	Immunohistochemistry: 1:50 -1:100 Dilutions listed as a recommendation. Optimal dilution should be determined by investigator.		
Storage:	Store at 4°C in the dark. Under these conditions, the antibody should have a shelf life of at least 12 months (provided they remain sterile). Do not freeze the antibody unless you want to store it for longer periods of time. Note, however, that each time an antibody preparation is frozen, about half its binding activity is lost.		

Application Notes

FOR RESEARCH USE ONLY

Neuromics' reagents are for in vitro and certain non-human in vivo experimental use only and not intended for use in any human clinical investigation, diagnosis, prognosis, or treatment. We disclaim all liability in connection with the use of the information contained herein or otherwise, and all such risks are assumed by the user.

v1