



Nop30

Datasheet

Catalog Number:	RA14126	Host:	Rabbit
Product Type:	Affinity purified antibody	Species Reactivity:	Human
Immunogen Sequence:	CPELQTLTRPGALRAPRRC Conjugated to KLH	Format:	Affinity Purified Liquid
Applications:	Western Blot 1:1000 Immunocytochemistry 1:500 Dilutions listed only as a recommendation. Optimal dilution should be determined by investigator.		
Storage:	Store frozen. Aliquot as undiluted antibody and immediately place at -20°C. Liquid may have become trapped in top of vial during shipping. Centrifugation of vial is recommended before opening. Stable for at least 6 months at -20°C. Repeated freeze/thaw cycles compromise the integrity of the antiserum.		

Application Notes

FOR RESEARCH USE ONLY

Neuromics' reagents are for in vitro and certain non-human in vivo experimental use only and not intended for use in any human clinical investigation, diagnosis, prognosis, or treatment. We disclaim all liability in connection with the use of the information contained herein or otherwise, and all such risks are assumed by the user.

v1