



ELISA Assay Diluents Standard Protocol

Catalog #: KF17347-51

Sizes: 25 ml, 100 ml, 500 ml, 1 L, 10 L

Storage: All of Neuromics' assay diluents include an antimicrobial agent to provide a useful shelf-life of 1 year when stored at 2-8C.

All proprietary assay diluents are compatible with most ELISA capture formats (however, high backgrounds may result when used with labeled generic anti-goat IgG signal generators). These diluents are designed to be compatible with Neuromics' plate coating and blocking reagents and should be matched with compatible sample matrices.

To use: Just pipette it directly onto the plate before you add the samples. We typically recommend using 50-100uL/well

FOR RESEARCH USE ONLY

v1-01/08

Neuromics' reagents are for in vitro and certain non-human in vivo experimental use only and not intended for use in any human clinical investigation, diagnosis, prognosis, or treatment. We disclaim all liability in connection with the use of the information contained herein or otherwise, and all such risks are assumed by the user.

phone 866-350-1500 • fax 612-677-3976 • www.neuromics.com