



Catalog Number:	RA14141	Host:	Rabbit
Product Type:	Whole Serum	Species	Rat
Immunogen Sequence:	SARPDAAKEEQQLRC	Reactivity:	
		Format:	100 ul. in PBS containing 0.05% sodium azide.
Applications:	Immunohistochemistry: 1:100-600 Dilutions listed as a recommendation. Optimal dilution should be determined by investigator.		
Storage:	Maintain at +2-8°C for 3 months or at -20°C for longer periods. Stable for 1 year. <i>Avoid repeated freeze-thaw cycles.</i>		

Application Notes

Immunohistochemistry:

Tested on cryo-preserved sections of the adult rat spinal cord (white matter) using anti-rabbit Cy3-conjugated secondary antibodies (red) and DAPI (blue) as a nuclear counter-stain.

General Reference:

Igor Jakovcevski and Nada Zecevic. Olig Transcription Factors Are Expressed in Oligodendrocyte and Neuronal Cells in Human Fetal CNS. *The Journal of Neuroscience*, November 2, 2005 25(44):10064 – 10073.