



Catalog Number: PR29005

Product Type: Purified Natural Protein

Source: Human Fibronectin Solution is manufactured exclusively from Human Plasma. Human Plasma (single donors) are tested for HBsAG-, HIV-, and HCV-antibodies and found negative.

Applications:

- ELISA
- Immunoblotting
- Cell culture substratum
- Stem cell research

Format: Sterile solution Buffered with 4M Urea and 1X TBS. pH=7.0-8.0. Concentration=1 mg/ml.

Storage: 2-8°C for up to 6 months. Aliquots may be stored in -80°C.

Avoid repeat freeze-thaw cycles.

Application Notes

****DO NOT AGITATE OR SWIRL**.** If entire amount of material is not to be used immediately, transfer into appropriate aliquots and store at -20°C for long term storage.

Advantages:

- Rapid attachment
- High plating and cloning efficiencies
- Good proliferation, high saturation density
- Lower requirements for serum and added growth factors
- Better response to physiologically occurring hormones, expression of differentiated functions
- Cells have longer life span, flattening and morphological changes, and improved plating consistency

Characteristics of Human Fibronectin:

- Human derived
- Sterile-filtered and cell-culture tested .
- Long-term storage when handled and stored properly under defined conditions.

FOR RESEARCH USE ONLY

NEUROMICS' REAGENTS ARE FOR IN VITRO AND CERTAIN NON-HUMAN IN VIVO EXPERIMENTAL USE ONLY AND NOT INTENDED FOR USE IN ANY HUMAN CLINICAL INVESTIGATION, DIAGNOSIS, PROGNOSIS, OR TREATMENT. THE ABOVE ANALYSES ARE MERELY TYPICAL GUIDES. THEY ARE NOT TO BE CONSTRUED AS BEING SPECIFICATIONS. ALL OF THE ABOVE INFORMATION IS, TO THE BEST OF OUR KNOWLEDGE, TRUE AND ACCURATE. HOWEVER, SINCE THE CONDITIONS OF USE ARE BEYOND OUR CONTROL, ALL RECOMMENDATIONS OR SUGGESTIONS ARE MADE WITHOUT GUARANTEE, EXPRESS OR IMPLIED, ON OUR PART. WE DISCLAIM ALL LIABILITY IN CONNECTION WITH THE USE OF THE INFORMATION CONTAINED HEREIN OR OTHERWISE, AND ALL SUCH RISKS ARE ASSUMED BY THE USER. WE FURTHER EXPRESSLY DISCLAIM ALL WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.-V1-11/2012

www.neuromics.com

Neuromics • 5325 West 74th Street, Suite 8 • Edina, MN 55439
phone 866-350-1500 • fax 612-677-3976 • e-mail pshuster@neuromics.com