



## Resistin


## Data Sheet

---

<b>Catalog Number:</b>	RA19059	<b>Host:</b>	Rabbit
<b>Product Type:</b>	Affinity Purified Anitbody	<b>Species</b>	Human
<b>Immunogen Sequence:</b>	Reacts with residues 94-102 [AGMDWTGAR] of the human Resistin protein.	<b>Reactivity:</b>	
		<b>Format:</b>	Liquid. 100 ug in 100 ul (1 mg/ml) in PBS containing 0.02% sodium azide
<b>Applications:</b>	<b>Western Blot:</b> 1:500-1,000 Dilutions listed as a recommendation. Optimal dilution should be determined by investigator.		
<b>Storage:</b>	Maintain at +2-8°C for 3 months or at -20°C for longer periods. Stable for 1 year. <i>Avoid repeated freeze-thaw cycles.</i>		
<b>Reference:</b>	Steppan C.M. et al. Proc Natl Acad Sci U S A 2001 Jan 16;98(2):502-6		

---

### Application Notes



FOR RESEARCH USE ONLY

Neuromics' reagents are for in vitro and certain non-human in vivo experimental use only and not intended for use in any human clinical investigation, diagnosis, prognosis, or treatment. We disclaim all liability in connection with the use of the information contained herein or otherwise, and all such risks are assumed by the user.

v1