

CLAC-P (NC2-2 region)

Data Sheet

Catalog Number: RA19008

Product Type: Rabbit IgG

Immunogen Sequence: KGEQGDQGPRMVFPK
(a.a. 155-169) of the NC2-2
region of human and mouse
CLAC-P

Format: Affinity Purified
liquid

Species Reactivity: Human

Applications: Western blotting

Dilutions: WB 1:500

Storage: Maintain at +2-8 degrees C for 3 months or at -20 degrees C for longer
periods. Stable for 1 year. Avoid repeated freeze-thaw cycles.
