



Catalog Number:	MO16102	Host:	Mouse
Product Type:	Protein G purified IgG ₁ Clone: 2H2.2	Species Reactivity:	Mouse
Immunogen Sequence:	Recombinant Murine LIF	Format:	Liquid 0.02M PBS, 0.25M NaCl pH 7.6 with 0.1% Sodium Azide
Applications:	Western blot: 1-10 µg/mL Dilutions listed as recommendation. Optimal dilution should be determined by the end user.		
Storage:	Maintain refrigerated at 2-8°C in undiluted aliquots for up to 12 months		

Application Notes

Reactivity to other LIF species is unknown.