

# NEUROMICS

## CERTIFICATE OF ANALYSIS

### Donor Donkey Serum

<b>Batch/Lot No:</b>	NA2041 (Formerly 206U20)	<b>Expiration Date:</b>	July 2025
<b>Catalog No:</b>	SER004	<b>Storage:</b>	≤ -10°C
<b>Country of Raw Material Origin:</b>	United States	<b>Appearance:</b>	Brown-Amber Liquid
<b>Country of Final Manufacturing:</b>	United States	<b>Filtration:</b>	Sterile 0.2 um membrane filtered
<b>Date of Manufacture:</b>	July 2020		

#### Chemical Parameters


Test	Methodology	Unit of Measure	Specification	Results
Hemoglobin	USP <90>	mg/dL	≤ 20	4.23
Osmolality	Osmometer	mOsm/kg H <sub>2</sub> O	260-340	277.0
pH	pH Meter		7.00-9.00	8.28
Total Protein	Chemistry Analyzer	g/dL	4.0-9.0	5.80

#### Microbiological Parameters

Test	Methodology	Unit of Measure	Specification	Results
Endotoxin	LAL	EU/mL	≤ 10	< 0.05
Mycoplasma	PCR		Not Detected	Not Detected
Sterility				
Bacteria	9CFR 113.26		No Growth	No Growth
Fungi	9CFR 113.26		No Growth	No Growth

**Test results reviewed and approved by\***

  
 \_\_\_\_\_  
 Approved by/Date      July 28, 2020

  
 \_\_\_\_\_  
 Verify by/Date      August 25, 2020

*\* Test results for this product may have been obtained from third-party certified testing facilities. Results are believed to be accurate but have not been verified by Neuromics. Neuromics relies on the integrity of third-party testing facilities and conducts periodic audits of those facilities. Such audits provided Neuromics with sufficient confidence in the facilities' level of compliance.*