

## Data Sheet

# GFP-Twin-Strep-tag<sup>®</sup> control protein

Cat. No.: 2-1007-005

Version: 2.2

Lot No.: 1007-

Revision Date: 25.02.2020

<b>Description</b>	GFP-Twin-Strep-tag <sup>®</sup> control protein for test purifications with Strep-Tactin <sup>®</sup> and Strep-Tactin <sup>®</sup> XT. The Twin-Strep-tag <sup>®</sup> is fused to the C-terminus of the green fluorescent protein (GFP).
<b>MW</b>	30 kDa
<b>Amount</b>	0.5 ml
<b>Concentration</b>	1 mg/ml
<b>Form</b>	dissolved in Strep-Tactin <sup>®</sup> /XT Wash Buffer; Buffer W
<b>Stability</b>	12 months after shipping
<b>Storage</b>	-25 to -15 °C
<b>Shipping</b>	room temperature
<b>Hazards</b>	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.



Go digital and help the environment. Please download all up-to-date manuals, protocols and other material from <http://www.iba-lifesciences.com>.

### For research use only

#### Important licensing information

This product is based on Twin-Strep-tag, Strep-Tactin technology covered by intellectual property (IP) rights and on completion of the sale IBA grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.iba-lifesciences.com/patents.html> or upon inquiry at [info@iba-lifesciences.com](mailto:info@iba-lifesciences.com) or at IBA GmbH, Rudolf-Wissell-Str. 28, 37079 Goettingen, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.

#### Trademark information

The owners of trademarks marked by "®" or "TM" are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.