

Data Sheet

Buffer XT-R

Strep-Tactin®XT Regeneration Buffer

Cat. No.: 2-1045-250

Lot No.: 1045-

Version: 1.0

Revision Date: 27.02.2020

Description	Buffer XT-R is a Strep-Tactin®XT Regeneration Buffer for competitive displacement of biotin from immobilized Strep-Tactin®XT.
Components	Buffer XT-R (Strep-Tactin®XT Regeneration Buffer, 250 ml): 3 M MgCl ₂
Stability	12 months after shipping
Storage	2–8°C
Shipping	room temperature
Hazards	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.
Working solution	The buffer is ready to use and should not be diluted.



Please download all up-to-date manuals, protocols and other material from
<http://www.iba-lifesciences.com>.

For research use only

Trademark information

The owners of trademarks marked by “®” or “TM” are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.