

Data Sheet

Strep-Tactin[®], lyophilized

Cat. No.: 2-1204-001, 2-1204-002, 2-1204-005,
2-1204-010, 2-1204-100

Version: 3.0
Revision Date: 06.07.2020

Lot No.: 1204 -

Description	Strep-Tactin [®] is a streptavidin variant for optimized Strep-tag [®] II binding.
Form	lyophilized from a 25 mg/ml solution in 10 mM potassium phosphate buffer pH 6.5
Extinction coefficient per subunit	$\epsilon_{280} = 35600 \text{ M}^{-1} \text{ cm}^{-1}$
MW per subunit	13342 Da
Stability (lyophilized)	6 months after shipping
Storage (lyophilized)	recommended: -25°C to -15 °C
Shipping	room temperature
Hazards	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.

Reconstitution	Dissolve with water
Storage (dissolved)	recommended: 2 – 8 °C, for long storage add 1 mM EDTA and/or 0.02 % NaN ₃ or pass the solution through a sterile filter



Please download all up-to-date manuals, protocols and other material from
<http://www.iba-lifesciences.com>.

For research use only

Important licensing information

This product is based on Strep-Tactin technology covered by intellectual property (IP) rights and on completion of the sale IBA grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.iba-lifesciences.com/patents.html> or upon inquiry at info@iba-lifesciences.com or at IBA GmbH, Rudolf-Wissell-Str. 28, 37079 Goettingen, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.

Trademark information

The owners of trademarks marked by “®” or “TM” are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.